Lampiran 1

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
TOSHIBA

1. OUTLINE

This is a capacitor discharge type mobile X-ray apparatus designed for clinical diagnostic use in hospitals and clinics.

This apparatus possesses the functions of automatic X-ray tube recharging and wave-tail cutoff, and enables radiographic control by the grid of a rotating anode. High safety, short-time radiographic operations can thus be performed on almost all sections of the patient's body. Also, the capacitor discharge system used in this apparatus eliminates the need for large power supply facilities and so allows direct use of the ordinary power supply that is provided to hospitals and clinics.

2. APPARATUS COMPOSITION


- (1) X-ray tube unit
- (2) High potential generator (high tension transformer and an X-ray control)
- (3) Mobile carriage (including an X-ray tube supporting mechanism)
- (4) Others (high potential cables, a hand switch)
- (5) Accessories

3. MAJOR OPERATIONAL CHARACTERISTICS

- (1) Maximum tube voltage : 100 kV
- (2) Tube current : DRX-66D ; 300 mA at 50 kV
150 mA at 100 kV
- (3) Effective focal spot size : 1,5mm
- (4) Target Angle : 15,5°
- (5) Anode : W, 80 mm Ø
- (6) High potential generation : Capacitor discharge method
- (7) Tube voltage adjustment range : 40 – 100 kV (continuously variable)
- (8) Capacitor charging time : Approx. 20 seconds from 0V to 100kV
- (9) Residual charge discharge method : Using the "DISCHARGE" push button on the control panel (Approx. 1 second from 100 kV to 0V)

ROTATING ANODE X-RAY
TUBE ASSEMBLY
CIRCLEX 1.5MG10AN-103 ✓
CIRCLEX 1.5MG10BN-10
(CONDENSER OPERATION)

03307

SHIMADZU CORPORATION
KYOTO JAPAN

1. SPECIFICATION

The CIRCLEX 1.5MG10AN-103 or the CIRCLEX 1.5MG10BN-10 is the grid controlled rotating anode X-ray tube assembly of miniature type. It is suitable for mass chest survey, the mobile X-ray unit and other low current radiography up to 300mA.

The X-ray tube assembly consists of the tube shield type RX-103 or RX10 with the X-ray tube CIRCLEX 1.5MG10AN or BN inserted.

Combined X-ray units are as follows.

X-ray tube assembly	X-ray unit
CIRCLEX 1.5MG10BN-10	SM-125L-20 (125kV)
CIRCLEX 1.5MG10AN-103 ✓	MC-100L-10 (100kV)

1.1 Technical Data

See Table 1.


Table 1

Max. operating voltage		kV _p	100(A), 125(B)
Effective focal spot size		mm	1.5
X-ray tube anode	Max. heat storage capacity	HU	60,000
	Max. cooling capacity	HU/min.	14,000
X-ray tube assembly	Max. heat storage capacity	HU	800,000
	Max. cooling capacity	HU/min.	12,000
Max. continu- ous load HU/sec.	X-ray tube	Anode stationary	240
		Anode rotating	260 ✓
	X-ray tube assembly		200
Max. radiographic rating		See radiographic rating chart.	
Filament characteristic		See filament characteristic chart.	
Emission characteristic		See emission characteristic chart.	
Cut-off grid voltage		V	-2000
Max. voltage between filament and grid		V	-3000
Target		Tungsten disc; target angle 18°; focal track diameter 40 mm	
Anode rotation		Anode should rotate anticlockwise, seen from cathode; revolution speed, over 2,700 rpm for 50Hz, over 3,200 rpm for 60 Hz,	
filtration		Approx. 1.5 mm Al	
X-ray field		Covered 35 cm X 35 cm at 65 cm from focus.	
X-ray protection		Less than 65 mR in an hour at one meter from focus except utilized beam.	
Min. ambient temperature		°C	-10
Weight		kg	Approx. 16
Stator		100V, 3.3A (starting) 40V, 1.4A (running) Starting time 0.8 sec. Split condenser 24 μF.	

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data indicated in this documentations refer to:

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PANTIX — rotating-anode tube

Standard-speed diagnostic rotating-anode tube with 2 focal spots of different size and rating. Preferred use of this rugged X-ray tube on generators of low and medium output.

The tungsten anode disk of 80 mm diameter allows the tube to be used for all routine examinations on fluoroscopic units, wall stands and planigraphic equipment as well as on Bucky examination tables.

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