LAMPIRAN

LAMPIRAN 1

PERHITUNGAN RALAT DENSITAS RADIOGRAF

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
LAMPIRAN 4

OPERATION MANUAL
FOR
BUCKY DEVICE
BF-10 (FOR BUCKY TABLE,
DIAGNOSTIC TABLE)
BF-11 (FOR BUCKY STAND)

PLEASE READ CAREFULLY
AND SAVE THIS MANUAL

SHIMADZU CORPORATION
KYOTO JAPAN

Table of contents

	Page
1. Application	1
2. Constitution	1
3. Specifications	1
3.1 Dimensions	1
3.2 Usable size of film	1
3.3 Grid	1
4. Operation	2
4.1 Connection of low voltage lead wire	2
4.2 Insertion of cassette	2
4.3 Replacement of Grid	2
5. Caution	3

1. Application

In combination with an X-ray high tension unit and an X-ray tube unit, this Bucky device BF-10 and BF-11 is used for Bucky radiography in vertical, horizontal and tilted positions. When used, this Bucky is generally incorporated in a vertical Bucky radiographic stand (in the case of BF-11), a Bucky table or a diagnostic table (in the case of BF-10).

2. Constitution

The Bucky BF-10 and BF-11 comprises the following parts:

- 2.1 Bucky main
- 2.2 Driving section
- 2.3 Switch section
- 2.4 Grid assembly
- 2.5 Tray

3. Specifications

The following are the specifications of Bucky BF-10 and BF-11.

- 3.1 Dimensions: 590 x 570 mm

- 3.2 Usable size of film:

Any size from 14 x 17 to 6½ x 8½ by lengthwise or crosswise, full exposure

- 3.3 Grid

Grid ratio R=8/1

Focal distance $f_0=100$ cm (in the case of Bucky table or diagnostic table Bucky BF-10)

$f_0=150$ cm (in the case of vertical Bucky radiographic stand Bucky BF-11)

Strips density: N=28 Lines/cm (BF-10)

N=34 Lines/cm (BF-11)

If required on a special order, other grids than described above are also available.

I. Operation (ref. Figs. 1 - 2)

4.1 Connection of low voltage lead wire

Connect terminals B₁, B₂, L₀, E and L₁₀₀ of the low voltage lead wire (4) to the corresponding terminals of the control of the X-ray unit.

4.2 Insertion of cassette

4.2.1 Draw out the cassette tray (2) by the grip (3)

4.2.2 Rotate the thumb piece (8) counterclockwise and the 2 cassette arresters (9) are loosened to make a gap suitable for insertion of a cassette.

Then, insert a cassette in the gap and while pressing the cassette firmly, lock it by rotating the thumb piece (8) clockwise.

Since these cassette arresters (9) are bilaterally adjustable, any size of cassette as mentioned in Item 3.2. above, when inserted in the gap between arresters, can be easily locked in the center position.

Inserting a cassette, lock it by thrusting it upward and put the spacer to down side of cassette in the case of a vertical Bucky radiographic stand, and lock the cassette after aligning the center line of the cassette with the (A) point of the V-indent at the center of cassette arresters (9) in other cases than a vertical Bucky radiographic stand. (ref. Fig. 2)

4.2.3 Return the cassette tray (2) to the original position.

4.2.4 In taking off the cassette tray (2), first draw it out as far as it stops and push the lock piece (10) and then take (2) off.

4.2.5 Insertion of cassette tray is possible from both sides of the bucky device.

4.3 Replacement of Grid

This Bucky device allows a grid to be replaced with easy operations.

The replacement of grid can be achieved in the following procedures.

4.3.1 Dimensions of replaceable grid

438 (lattice-directionally) x 479 mm

4.3.2 Turn the line switch on.

4.3.3 Remove the cassette tray. (See 4.2.4.)


Pull the knob (6) and the rails for cassette tray, on which the grid is mounted, are inclined. This allows the motor to run and the grid to move in your direction.

- 4.3.4 When the lever (11) is raised up the grid comes out by the spring action.
Remove the grid and replace it to another.
Care should be taken not to reverse the grid.
- 4.3.5 Insert the replaced grid to the rail for grid, and make sure that the grid is across the stopper (12).
- 4.3.6 As soon as the grid is raised up, the grid is reset in the original position.
- 4.3.7 Set the cassette tray in a previous manner.

5. Caution


In advance to assembling the table top, be sure to remove the screw that is fixing the grid of Bucky unit.

- | | | |
|---------------------------|-----------------------|-----------------|
| (1) Bucky frame | (2) Cassette tray | (3) Grip |
| (4) Low voltage lead wire | (6) Knob | (7) Grid holder |
| (8) Thumb piece | (9) Cassette arrester | (10) Lock piece |
| (11) Lever | (12) Stopper | |

Fig. 1 Whole view

(for use on vertical radiographic stand) (for use on Bucky table or diagnostic table)

Fig. 2 Cassette rock

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