## DAFTAR PUSTAKA

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● Indicates analyzer(s) on which cobas c pack can be used

COBAS INTEGRA 400/400 plus	COBAS INTEGRA 800
●	●

**System information**  
COBAS INTEGRA Cholesterol Gen.2 (CHOL2)  
Test CHOL2, test ID 0-586

**Intended use**  
In vitro test for the quantitative determination of total cholesterol in serum and plasma on COBAS INTEGRA systems.

**Summary**<sup>1,2,3,4,5,6,7,8,9,10,11,12</sup>  
Cholesterol is a steroid with a secondary hydroxyl group in the C3 position. It is synthesized in many types of tissue, but particularly in the liver and intestinal wall. Approximately three quarters of cholesterol are newly synthesized and a quarter originates from dietary intake. Cholesterol assays are used for screening for atherosclerotic risk and in the diagnosis and treatment of disorders involving elevated cholesterol levels as well as lipid and lipoprotein metabolic disorders. Cholesterol analysis was first reported by Liebermann in 1885 followed by Burchard in 1889. In the Liebermann-Burchard reaction, cholesterol forms a blue-green dye from polymeric unsaturated carbohydrates in an acetic acid/acetic anhydride/concentrated sulfuric acid medium. The Abell and Kendall method is specific for cholesterol, but is technically complex and requires the use of corrosive reagents. In 1974, Roeschlau and Allain described the first fully enzymatic method. This method is based on the determination of  $\Delta^4$ -cholestenone after enzymatic cleavage of the cholesterol ester by cholesterol esterase, conversion of cholesterol by cholesterol oxidase,

and subsequent measurement by the Trinder reaction of the hydrogen peroxide formed. Optimization of ester cleavage (> 99.5 %) allows standardization using primary and secondary standards and a direct comparison with the CDC and NIST reference methods. Nonfasting sample results may be slightly lower than fasting results. The Roche cholesterol assay meets the 1992 National Institutes of Health (NIH) goal of less than or equal to 3 % for both precision and bias.

**Test principle**  
Enzymatic, colorimetric method  
Cholesterol esters are cleaved by the action of cholesterol esterase to yield free cholesterol and fatty acids. Cholesterol oxidase then catalyzes the oxidation of cholesterol to cholest-4-en-3-one and hydrogen peroxide. In the presence of peroxidase, the hydrogen peroxide formed effects the oxidative coupling of phenol and 4-aminoantipyrine to form a red quinone-imine dye.

$$\text{Cholesterol esters} + \text{H}_2\text{O} \xrightarrow{\text{CE}} \text{cholesterol} + \text{RCOOH}$$

$$\text{Cholesterol} + \text{O}_2 \xrightarrow{\text{CHOD}} \text{cholest-4-en-3-one} + \text{H}_2\text{O}_2$$

$$2 \text{H}_2\text{O}_2 + 4\text{-AAP} + \text{phenol} \xrightarrow{\text{POD}} \text{quinone-imine dye} + 4 \text{H}_2\text{O}$$

The color intensity of the dye formed is directly proportional to the cholesterol concentration. It is determined by measuring the increase in absorbance at 512 nm.

2009-05, V 4 EN

1 / 4

**CHOL2**

**Reagents - working solutions**

Components	Concentrations		
	R	Test	
PIPES <sup>a</sup> buffer	225	75	mmol/L
Mg <sup>2+</sup>	10	3.3	mmol/L
Sodium cholate	0.6	0.2	mmol/L
4-Aminoantipyrine	≥ 0.45	≥ 0.15	mmol/L
Phenol	≥ 12.6	≥ 4.2	mmol/L
Fatty alcohol polyglycol ether	3	1	%
CE (Pseudomonas spec.)	≥ 1.5	≥ 0.5	kU/L
		(≥ 8.3	μkat/L)
CHOD (E. coli)	≥ 0.45	≥ 0.15	kU/L
		(≥ 2.5	μkat/L)
POD (horseradish)	≥ 0.75	≥ 0.25	kU/L
		(≥ 4.1	μkat/L)
pH	6.8	6.8	

a) PIPES = Piperazine-1,4-bis(2-ethanesulfonic acid)

The reagent contains nonreactive preservative and stabilizer.

**Precautions and warnings**

Pay attention to all precautions and warnings listed in this Method Manual, Chapter 1, Introduction.

**Reagent handling**

Ready for use.

**Storage and stability**

Shelf life at 2 to 8 °C See expiration date on cobas c pack label.

COBAS INTEGRA 400/400 plus systems

On-board in use at 10 to 15 °C 8 weeks

COBAS INTEGRA 800 systems

On-board in use at 8 °C 8 weeks

**Specimen collection and preparation**

For specimen collection and preparation, only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable. Serum.

Plasma: Li-heparin or K<sub>3</sub>-EDTA plasma. (Use of EDTA-plasma leads to slightly lower values.)

Do not use citrate, oxalate, or fluoride.<sup>13</sup>

Fasting and nonfasting samples can be used.<sup>11</sup>

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Stability:<sup>14,15</sup> 7 days at 15-25 °C  
7 days at 2-8 °C  
3 months at (-15)-(-25) °C

Centrifuge samples containing precipitates before performing the assay.

**Materials provided**

See "Reagents - working solutions" section for reagents.

**Assay**

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

**Application for serum and plasma****COBAS INTEGRA 400/400 plus test definition**

Measuring mode	Absorbance
Abs. calculation mode	Endpoint
Reaction mode	R-S
Reaction direction	Increase
Wavelength A/B	512/659 nm
Calc. first/last	17/69
Unit	mmol/L

**Pipetting parameters**

		Diluent (H <sub>2</sub> O)
R	47 μL	70 μL
Sample	2 μL	23 μL
Total volume	142 μL	

**COBAS INTEGRA 800 test definition**

Measuring mode	Absorbance
Abs. calculation mode	Endpoint
Reaction mode	R-S
Reaction direction	Increase
Wavelength A/B	512/659 nm
Calc. first/last	17/98
Unit	mmol/L

**Pipetting parameters**

		Diluent (H <sub>2</sub> O)
R	47 μL	73 μL
Sample	2 μL	20 μL
Total volume	142 μL	

**Calibration**

Calibrator	Calibrator f.a.s. Use deionized water as zero calibrator.
Calibration mode	Linear regression
Calibration replicate	Duplicate recommended
Calibration interval	Each lot and as required following quality control procedures

Traceability: This method has been standardized by ID-MS<sup>b</sup> and also according to Abell-Kendall.

This complies with the requirements of the National Institute of Standards and Technology (NIST).

b) Isotope dilution - mass spectrometry

**Quality control**

Reference range	Precinorm U, Precinorm U plus or Precinorm L
Pathological range	Precipath U, Precipath U plus or Precipath L
Control interval	24 hours recommended
Control sequence	User defined
Control after calibration	Recommended

For quality control, use control materials as listed in the Order information section. Other suitable control material can be used in addition.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits.

Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Follow the applicable government regulations and local guidelines for quality control.

**Calculation**

COBAS INTEGRA analyzers automatically calculate the analyte concentration of each sample. For more details, please refer to Data Analysis in the Online Help (COBAS INTEGRA 400/400 plus/800 analyzers).

Conversion factor:  $\text{mmol/L} \times 38.66 = \text{mg/dL}$

**Limitations - interference<sup>16</sup>**

Criterion: Recovery within  $\pm 10\%$  of initial value.

*Serum, plasma*

Icterus	No significant interference up to an I index of 16 for conjugated bilirubin and 11 for unconjugated bilirubin (approximate conjugated bilirubin concentration: 274 $\mu\text{mol/L}$ or 16 mg/dL, and approximate unconjugated bilirubin concentration: 188 $\mu\text{mol/L}$ or 11 mg/dL). <sup>c</sup>
Hemolysis	No significant interference up to an H index of 810 (approximate hemoglobin concentration: 503 $\mu\text{mol/L}$ or 810 mg/dL). <sup>c</sup>
Lipemia	No significant interference. <sup>c</sup>
Drugs	No interference was found at therapeutic concentrations using common drug panels. <sup>17,18</sup>
Other	In very rare cases gammopathy, in particular type IgM (Waldenström's macroglobulinemia), may cause unreliable results.

c) measured at cholesterol levels up to 5.28 mmol/L

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

**ACTION REQUIRED**

**Special wash programming:** The use of special wash steps is mandatory when certain test combinations are run together on COBAS INTEGRA analyzers. Refer to the Method Manual, Introduction, Extra Wash Cycles for further instructions. Where required, special wash/carry-over evasion programming must be implemented prior to reporting results with this test.

**Limits and ranges****Measuring range**

0.1-20.7 mmol/L (3.87-800 mg/dL)

Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:10 dilution. Results from samples diluted by the rerun function are automatically multiplied by a factor of 10.

**Lower limits of measurement**

Lower detection limit of the test:

0.1 mmol/L (3.87 mg/dL)

The detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying three standard deviations above that of a zero sample (zero sample + 3 SD, within-run precision, n = 21).

**Expected values**

Clinical interpretation according to the recommendations of the European Atherosclerosis Society:<sup>19</sup>

	mmol/L	mg/dL	Lipid metabolic disorder
Cholesterol	< 5.2	< 200	No
Triglycerides	< 2.3	< 200	
Cholesterol	5.2-7.8	200-300	Yes, if HDL-cholesterol < 0.9 mmol/L (< 35 mg/dL)
Cholesterol	> 7.8	> 300	Yes
Triglycerides	> 2.3	> 200	

Recommendations of the NCEP Adult Treatment Panel for the following risk-cutoff thresholds for the US American population:<sup>20</sup>

Desirable cholesterol level	< 5.2 mmol/L	(< 201 mg/dL)
Borderline high cholesterol	5.2-6.2 mmol/L	(200-240 mg/dL)
High cholesterol	$\geq 6.2$ mmol/L	( $\geq 240$ mg/dL)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

**Specific performance data**

Representative performance data on COBAS INTEGRA analyzers are given below. Results obtained in individual laboratories may differ.

**Precision**

Precision was determined using human samples and controls in an internal protocol. Repeatability<sup>d</sup> (n = 21), intermediate precision<sup>e</sup> (1 aliquots per run, 1 run per day, 21 days). The following results were obtained:

	Level 1	Level 2
Mean	2.74 mmol/L (106 mg/dL)	6.20 mmol/L (240 mg/dL)
CV repeatability <sup>d</sup>	0.51 %	0.81 %
Mean	2.61 mmol/L (101 mg/dL)	5.96 mmol/L (230 mg/dL)
CV intermediate precision <sup>e</sup>	1.9 %	1.4 %

d) repeatability = within-run precision

e) intermediate precision = total precision/between-run precision/between-day precision

**Method comparison**

Cholesterol values for human serum and plasma samples obtained on a COBAS INTEGRA 700 analyzer with the COBAS INTEGRA Cholesterol Gen.2 reagent (CHOL2) were compared to those determined by ID-MS.

ID-MS	Sample size (n) = 50
Passing/Bablok <sup>21</sup>	Linear regression
$y = 0.99x + 0.04$ mmol/L	$y = 0.98x + 0.09$ mmol/L
$r = 0.971$	$r = 0.999$
SD (md 95) = 0.115	$Sy.x = 0.058$

The sample concentrations were between 1.51 and 10.94 mmol/L (58.4 and 423 mg/dL).

**References**

- Greiling H, Gressner AM, eds. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd ed. Stuttgart/New York: Schattauer, 1995.
- Liebermann C. Ber Dtsch Chem Ges 1885;18:1803.
- Burchard H. Beiträge zur Kenntnis der Cholesterine. Dissertation, Rostock 1889.
- Abell L et al. Standard Methods in Clinical Chemistry 1958;26:2.
- Allain CC et al. Clin Chem 1974;20:470.
- Roeschlau P et al. Z Klin Chem Klin Biochem 1974;12:226.
- Trinder P. Ann Clin Biochem 1969;6:24.
- Siedel J, Hägele EO, Ziegenhorn J et al. Clin Chem 1983;29:1075.
- Wiebe DA, Bernert JT. Clin Chem 1984;30:352.
- Cohn JS, McNamara JR, Schaefer EJ. Lipoprotein Cholesterol Concentrations in the Plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459.
- Pisani T, Gebiski CP, Leary ET et al. Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995;119:1127.

## LAMPIRAN 7. Hasil uji statistik

### 1. Jumlah sampel tiap kelompok

Case Processing Summary						
Kelompok	Cases					
	Valid		Missing		Total	
	N	Percent	N	Percent	N	Percent
Kadar Kolesterol Kontrol	5	100.0%	0	.0%	5	100.0%
Minyak	5	100.0%	0	.0%	5	100.0%
Ubi	5	100.0%	0	.0%	5	100.0%
Minyak + Ubi	5	100.0%	0	.0%	5	100.0%

### 2. Analisis deskriptif pengaruh pemberian ubi ungu terhadap kadar kolesterol total serum

Descriptives			Statistic	Std. Error
Kelompok				
Kadar Kolesterol	Kontrol	Mean	57.20	5.962
		95% Confidence Interval for Mean		
		Lower Bound	40.65	
		Upper Bound	73.75	
		5% Trimmed Mean	56.61	
		Median	54.00	
		Variance	177.700	
		Std. Deviation	13.330	
		Minimum	45	
		Maximum	80	
		Range	35	
		Interquartile Range	19	

	Skewness		1.733	.913
	Kurtosis		3.608	2.000
Minyak	Mean		59.60	4.020
	95% Confidence Interval for Mean	Lower Bound	48.44	
		Upper Bound	70.76	
	5% Trimmed Mean		59.33	
	Median		56.00	
	Variance		80.800	
	Std. Deviation		8.989	
	Minimum		52	
	Maximum		72	
	Range		20	
	Interquartile Range		17	
	Skewness		.714	.913
	Kurtosis		-1.842	2.000
Ubi	Mean		55.60	4.665
	95% Confidence Interval for Mean	Lower Bound	42.65	
		Upper Bound	68.55	
	5% Trimmed Mean		55.33	
	Median		52.00	
	Variance		108.800	
	Std. Deviation		10.431	
	Minimum		45	
	Maximum		71	
	Range		26	
	Interquartile Range		19	
	Skewness		.839	.913
	Kurtosis		-.420	2.000

Minyak + Ubi	Mean		54.80	4.067
	95% Confidence Interval for Mean	Lower Bound	43.51	
		Upper Bound	66.09	
	5% Trimmed Mean		54.39	
	Median		51.00	
	Variance		82.700	
	Std. Deviation		9.094	
	Minimum		47	
	Maximum		70	
	Range		23	
	Interquartile Range		15	
	Skewness		1.592	.913
	Kurtosis		2.568	2.000

3. Hasil uji normalitas data pengaruh pemberian ubi ungu terhadap kadar kolesterol total serum

**Tests of Normality**

Kelompok	Kolmogorov-Smirnov <sup>a</sup>			Shapiro-Wilk		
	Statistic	df	Sig.	Statistic	df	Sig.
Kadar Kolesterol Kontrol	.366	5	.028	.809	5	.097
Minyak	.256	5	.200 <sup>*</sup>	.860	5	.228
Ubi	.235	5	.200 <sup>*</sup>	.935	5	.634
Minyak + Ubi	.262	5	.200 <sup>*</sup>	.847	5	.184

a. Lilliefors Significance Correction

\*. This is a lower bound of the true significance.

## 4. Hasil uji homogenitas varians

Test of Homogeneity of Variance

		Levene Statistic	df1	df2	Sig.
Kadar Kolesterol	Based on Mean	.178	3	16	.910
	Based on Median	.054	3	16	.983
	Based on Median and with adjusted df	.054	3	13.299	.983
	Based on trimmed mean	.136	3	16	.937

5. Hasil uji parametrik *One Way Anova*

ANOVA

Kadar Kolesterol

	Sum of Squares	df	Mean Square	F	Sig.
Between Groups	67.200	3	22.400	.199	.895
Within Groups	1800.000	16	112.500		
Total	1867.200	19			

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